TIGER
D E N T A L

Innovations in orthodontics

TOPJET *clix*

Gebrauchsanweisung
User Manual

TopJet clix Distalizer

250, 360, 250+8, 360+8

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
TopJet 250 = 40 g/Raste, TopJet 360 = 60 g/Raste

English

Dear doctor, we are honored that you have chosen this innovative new product from TigerDental.

To ensure the safe and effective use of this product, you must carefully read and adhere to this user manual.

The user manual covers all necessary safety instructions but cannot describe all possible circumstances. For this reason, you may find further user notices online at tigerdental.com under the “Know-how” link. Our customer service is available for any additional questions.

The continual development of our products indicates our company’s commitment to innovation. Accordingly, after several clinical uses, we recommend you carefully review the accompanying user manual (this can be found online under the “Downloads” link at tigerdental.com).

1. Manufacturer

TigerDental, Belruptstrasse 59, 6900 Bregenz, Austria

2. Product Description

TopJet clix allows for reliable and incremental molar distalization through the use of a programmable telescoping arm. The operator benefits from the appliance’s optimized handling characteristics and patients benefit from its maximal comfort. The TopJet clix appliance offers appealing esthetics and enables compliance free one or two-sided distalization of the desired teeth of up to 0.8 mm per month. Due to its compact design, the installation of TopJet clix is amazingly simple and can be done in small and narrow palates. The TopJet has already been successfully used in patients 8 years and older. The appliance can be installed in just a few minutes. Add this feature to near invisibility and high patient comfort, and it is easy to see why this is such a great device. During the entire installation process, a thread secures the TopJet Distalizer to its packaging outside the mouth. This safety feature minimizes the risk of aspiration. The appliance can only be activated after the tread is cut. The appliance’s hinge-like connection to the palatal arch, called the T-connector, prevents distal tipping of the molars. Use of a transpalatal arch (TPA) ensures a lack molar rotation and the TopJet’s friction-free power output. Operators can choose from either a

250 cN (black coding) or 360 cN (red coding) version (see section “Selecting the TopJet clix”).

The TopJet clix can be special ordered for lengths in excess of 8 mm. At the maximum spring activation, the smallest distance between the micro implant and the TPA is 14 mm while the largest distance is 28 mm. This allows for a maximum distalization of 14 mm. The TopJet can be re-activated in 1.8 mm increments.

Selection of TopJet clix

Black Coding: TopJet clix 250 (250 cN) Is mainly used in children to distalize the upper first molars only.

Red Coding: TopJet clix 360 (360 cN) For molars and premolar distalization. Can be deactivated to exert less force.

3. Field of Application

Indications

1. Crowding in the maxillary anterior with sufficient space in the posterior for distalization (consider removal of third molars!)
2. Protruded maxillary anterior teeth
3. Lack of space in the maxillary buccal segments (i.e. peg laterals, missing/narrow primary teeth or impacted premolars).
4. Class II malocclusion with mesially migrated maxillary molars (also in skeletal class III/dental class II combination cases). The TopJet clix can be used without problems on one side or both sides.

Contraindications

1. Periodontally compromised molars when using high force values
2. Short maxilla (distally directed second molars)

4. Installation of TopJet clix

System Prerequisites

- 16/26 bands with windowed Goshgarian lingual sheaths
- TopJet transpalatal arch (prefabricated or using bending template)
- Dual Top micro implant (JetScrew) measuring: 2.0 mm x 10/12/14 mm, 2.5 mm x 12/14/16 mm and placed in the M4 position
- Off-angle KaVo ENDOstraight handpiece with a 25 mm hex driver attachment
- A hooked explorer for activating and deactivating the TopJet
- Ligature cutter to cut the security thread
- TopJet-Pliers: to close the T- connector
- Tweezers: holding the TopJet clix during installation

Insertion the new TopJet clix

The C-clip is secured to the head of the micro implant with a light force. Afterward, tweezers are used to bring TopJet clix into contact with the TPA. Using the T-connector pliers (modified Tweed pliers) to close the flaps of the T-connector allowing for some play. After locking the TopJet clix between the micro implant and the TPA apply a light-curing resin (Triad® Gel) to each of these connections. This creates a rotationally stable connection between screw and TopJet clix (Fig. 1–5) and increases patient comfort.

Activation of TopJet clix

Now activate TopJet clix by cutting the security thread between the ball and the spring telescope. The patient will now feel a slight pressure now on the molars (Fig. 6).

Easy Reactivation and Deactivation

Activation as well as deactivation are possible in seconds with TopJet clix using a hooked explorer or probe. During molar distalization, spring tension can be increased by hooking a probe behind the distal end of the telescope tube and pulling mesially. The spring will lock into place audibly. Too much activation? To deactivate the telescope simply disengage the spring lock with the same hooked explorer. When disengaging the spring lock, both tubes will slide distally, simultaneously decreasing the force on the molars. Fine adjustment in this fashion is extremely easy and can be selected at increments of 1.8 mm.

1 Insertion

2 Completed insertion

3 Slide the clip on the screw

4

5 Extend the T-connector

6 Position the T-connector


7 Close the T-connector

8 Coat the T-connector

9 Coat the screw head

10 Activation: Remove security thread

11 Activation/reactivation

12 Deactivation

13 Traditional Removal: Press out the TPA from the distal side

14 In Case of Emergency: Cutting through the TPA

15 Remove Plastic-coat

16 Slide the c-clip off the screw

17 C-clip is now detached

18 Detaching the T-connector

19 Swivel the TopJet upwards to loosen/detach the T-connector

20 Remove the TopJet

5. Removal of TopJet clix

Traditional Removal of the TPA

A TPA inserted as per the previously described procedure will adhere to the Goshgarian lingual sheath with great strength. The easiest way remove the TPA is to simultaneously apply force to the to distally protruding portion of the TPA and the lingual sheath. Alternatively, the 2-plier method is feasible. For this purpose, a ligature cutter and a WEINGART plier are used. The front ends of each of these instruments are placed in close proximity so that they can pivot when the handles of the instruments are moved closer together. The resulting lever action removes the TPA from the Goshgarian lingual sheath. In Case of Emergency: Cutting Through the TPA

In the event that a traditional removal of the TPA is not possible (i.e. pressure from the TPA has ulcerated the palate), it is recommended to use the small diamond disc in an off-angle handpiece with water cooling to carefully sever the TPA. During this process ensure that the TPA is secured with mosquito pliers, needle-nosed pliers or with floss. Removal of the C-clip from the Micro Implant Prior to removal of the C-clip, all Triad® resin should be removed from the micro implant head with a ligature cutter. The C-clip at the front end of the TopJet can usually be removed easily with WEINGART pliers. Holding one arm of the plier to the micro implant head, push the protruding end of the C-clip toward the neck of the micro implant. Any remaining Traid® gel adhesive residues will be broken in this process. Attention! Do not remove the C-clip directly from the micro implant by “pulling” on the spring pin. This risks the screw loosening! Unlocking the T-connector and removing the TopJet clix. After removing C-clip from the micro implant, grasp the spring pin with WEINGART pliers. Use a 45° rotary movement to open the flaps of the T-connectors and to break the adhesive connection to the TPA. Now the TopJet clix can be easily removed.

6. Safety Instructions

This operating and installation manual contains tips and recommendations, but does not release you from personal responsibility! The use of the TopJet clix is allowed only by orthodontists and dentists accord-

ing to Medical Device Directive (MDD). The user must comply with the provisions of the Medical Devices Directive and any other applicable laws and regulations. Patients with an allergy to nickel and/or chromium should not use this product. The doctor should be consulted immediately regarding any irritation to suspend treatment.

7. Notices for Single-use Products

The TopJet clix is nearly ready to use out of the box, but is not sterile and should not be used more than once. The reprocessing of a once used TopJet clix (recycling) and its reapplication to the patient is strictly prohibited.


8. References of Quality


Tiger Dental assures the user a flawless quality of its products. The content of this manual is based on experience. The user is responsible for the processing of the products themselves. TigerDental is not liable for erroneous results or processing by the user that occur in the absence of influence of Tiger Dental. Although the TopJet clix is a mature product, improvements and changes can be made, which may not be documented in the present manual. Similarly, we are not liable for printing errors.

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